

DatamedWL™ Compatibility List

Updated April, 2025 – DatamedWL™ v3.2.0

Devices

Baxter / Hillrom (Mortara®, Welch Allyn®)

Devices: ELI® series cardiographs and Welch Allyn® Diagnostic Cardiology Suite (DCS) (formerly Connex Cardio)

Transmission Method(s): LAN/WLAN

Notes: The ELI Link utility (available from Mortara) is required, and the cardiograph must have the DICOM® option. ELI Link must be configured to retrieve orders from DatamedWL via DICOM Modality Worklist, and to send ECGs to Datamed via DICOM ECG Waveform or Mortara XML files. For Server 2016/2019 the latest version of ELI Link is needed.

Website: www.hillrom.com

Devices: Welch Allyn® CP150 with DICOM option

Transmission Method(s): (W)LAN

Notes: CP150 with optional DICOM capability can communicate bidirectionally with worklist download and ECG waveform upload. Optional wireless radio allows for wireless DICOM communication.

Website: www.hillrom.com

Devices: Burdick® Atria 3100/6100 and 8300/8500 cardiographs

Transmission Method(s): LAN/WLAN

Notes: The cart must have an Ethernet/wireless adapter for communications. Serial/modem connections are not supported. Atria 3000 has no order capability and is only compatible with DatamedFT, not DatamedWL.

Devices: Eclipse™ Premier cardiographs

Transmission Method(s): LAN/WLAN

Notes: The cart must have the Eclipse Wireless Client Bridge (WCB) option. Serial/modem connections are not supported. **IMPORTANT: Order Number is not included in the EKG that is sent from these devices. Some systems (e.g. Epiphany Cardio Server® and GE® MUSE®) can associate and close out the order without the Order Number, but other systems (e.g. McKesson Cardiology™) require Order Number to close out the order. Keep this in mind when planning the order workflow.**

DICOM®

Devices:

Any device that supports DICOM Modality Worklist

Notes:

Order query will use DICOM MWL. ECGs must be sent as DICOM waveform objects via a C-STORE transaction.

Edan Instruments, Inc.

Devices: SE-12 Express, SE-1200 Express, SE-1201, SE-1202, SE-15, SE-1515, SE-16, SE-18, and SL-18 cardiographs

Transmission Method(s): Network transmission (LAN+FTP)

Notes:

The cart should be set to SCP (.scp) output. An FTP server such as the Windows Server® FTP server or the open source program FileZilla is required.
IMPORTANT: Order Number is not included in the EKG that is sent from these devices. Some systems (e.g. Baxter Cardio Server® and GE® MUSE®) can associate and close out the order without the Order Number, but other systems (e.g. Change Healthcare Cardiology™) require Order Number to close out the order. Keep this in mind when planning the order workflow.

Website:

www.edan.com

Fukuda Denshi

Devices: FX-8400 cardiographs

Transmission Method(s): LAN/WLAN

Notes:

Communicates using DICOM MWL configuration.

Website:

www.fukuda.com

GE® Healthcare

Devices: MAC® VU360, MAC® 5, MAC® 7, and MAC® 2000 cardiographs

Transmission Method(s): LAN/WLAN (DCP/Web)

Website:

www.gehealthcare.com

Devices: MAC® 800, MAC® 2000, MAC® 3500, MAC® 5000, MAC® 5500, and MAC® 5500HD cardiographs

Transmission Method(s): LAN/WLAN (CSI)

Notes:

The Order Manager must be enabled on the cardiograph. The v9 firmware can only query using patient id where v10 can use visit id, but most systems can't handle the Visit ID query.

Website:

www.gehealthcare.com

Nihon Kohden® Corporation

Devices: Cardiofax® S (ECG-1250, ECG-2250, ECG-3250)

Cardiofax® M (ECG-1350, ECG-2350, ECG-3350)

Cardiofax® V (ECG-1500, ECG-1550)

Cardiofax® G (ECG-2450, ECG-2550)

Transmission Method(s): Network transmission (LAN+FTP)

Notes: Datamed LLC will supply NK's receiver utility. NK has updated the utility which is not compatible with Server 2016/2019. An FTP server is also required – the open source program FileZilla has been tested.

Website: www.nihonkohden.com

Philips® Healthcare

Devices: PageWriter® Touch, PageWriter® TC cardiographs

Transmission Method(s): LAN/WLAN

Website: www.usa.philips.com

SCHILLER AG

Devices: CARDIOVIT® AT-10 plus, AT-102, AT-102 Plus, AT-170, FT-1, AT-102 G2, MS-2007, MS-2010, and MS-2015 cardiographs

Transmission Method(s): LAN/WLAN

Notes: These devices have worklist available as an option from SCHILLER. The worklist option in AT-102 also requires the built-in SCM option.

Website: www.schiller.ch

Devices: CARDIOVIT® CS-104, CS-200, CS-200 Office, CS-200 Touch, and CS-200 Excellence cardiographs

Transmission Method(s): LAN/WLAN

Notes: These devices have worklist built-in.

Website: www.schiller.ch

Systems

AGFA HealthCare Corporation

System Name: AGFA HealthCare Enterprise Imaging™
Notes: Supports DICOM MWL for order query. ECGs are sent as DICOM waveform objects via a C-STORE transaction.
Website: global.agfahealthcare.com

Baxter / Hillrom (Epiphany, Mortara®)

System Name: Baxter Cardio Server
Notes: Supports DICOM MWL for order query. DatamedFT™ will be configured to write the ESI format ECG files to the Cardio Server input folder.
Website: Baxter Cardio Server
System Name: Pyramis®
Notes: For this system DatamedFT™ will be configured to write the converted ECG files to the system server.

Carestream Health, Inc.

System Name: Carestream® Cardiology PACS
Notes: Supports DICOM MWL for order query. ECGs are sent as DICOM waveform objects via a C-STORE transaction. In the Storage/Management System Settings set System Type to ScImage – PicomEnterprise.
Website: www.carestream.com

Cerner Corporation (Oracle)

System Name: Powerchart® ECG
Notes: Supports DICOM MWL for order query. ECGs are sent as DICOM waveform objects via a C-STORE transaction.
Website: www.oracle.com/health

Change Healthcare (Optum)

System Name: Change Healthcare Cardiology ECG Management™
Notes: Supports DICOM MWL for order query. DatamedFT™ will be configured to write the MCK format ECG files to the system server.
Website: business.optum.com
System Name: McKesson Cardiology ECG Management™
Notes: For this system DatamedFT™ will be configured to write the MCK format ECG files to the system server.

DICOM®

System Name:

Any system that supports DICOM Modality Worklist

Notes:

Order query will use DICOM MWL. ECGs can be sent as DICOM waveform objects via a C-STORE transaction or as a file in any format supported by DatamedFT™.

Esaote®

System Name:

Suitestensa

Notes:

For this system DatamedFT™ will be configured to write the ESA format ECG files to the system server. Alternatively the ECGs can be sent as DICOM waveform objects via a C-STORE transaction.

Website:

www2.esaote.com

Fujifilm Medical Systems

System Name:

VidiStar® PACS and Online Reporting Software

Notes:

Supports DICOM MWL for order query. ECGs are sent as DICOM waveform objects via a C-STORE transaction.

Website:

www.fujimed.com

System Name:

Synapse® Cardiovascular

Notes:

Supports DICOM MWL for order query. DatamedFT™ will be configured to write the ESI format ECG files to the Synapse ECG input folder. Alternatively ECGs can be sent as DICOM waveform objects via a C-STORE transaction.

Website:

www.fujimed.com

GE® Healthcare

System Name:

MUSE®

Notes:

The MUSEAPI3™ option is required and can be installed on v7.1.1 or later.

Website:

www.gehealthcare.com

INFINITT Healthcare Co

System Name:

INFINITT Cardiology PACS

Notes:

Supports DICOM MWL for order query. ECGs are sent as DICOM waveform objects via a C-STORE transaction.

Website:

www.infinitt.com

Intelerad (formerly LUMEDX Corporation)

System Name: HealthView ECG Manager

Notes: Supports DICOM MWL for order query. DatamedFT™ will be configured to write the XML format ECG files to the server.

Website: www.intelerad.com

Medimatic

System Name: ComPACS™

Notes: Supports DICOM MWL for order query. ECGs are sent as DICOM waveform objects via a C-STORE transaction.

Website: www.medimatic.com

Neagen Oy

System Name: Neagen PACS

Notes: Supports DICOM MWL for order query. ECGs are sent as DICOM waveform objects via a C-STORE transaction.

Website: www.neagen.com

Novarad® Corporation

System Name: NovaCardio® ECG

Notes: Supports DICOM MWL for order query. ECGs are sent as DICOM waveform objects via a C-STORE transaction.

Website: www.novarad.net

Philips® Healthcare

System Name: IntelliSpace® ECG

Notes: The iECG system must have IntelliBridge Enterprise (IBE), which DatamedWL™ will communicate with directly. DatamedFT™ will be configured to write the XML format ECG files to the IntelliSpace® ECG server.

Website: www.usa.philips.com

System Name: TraceMasterVue™

Notes: DatamedWL™ communicates directly with Philips OrderVue™.

Website: www.usa.philips.com

ScImage®

System Name:

PicomEnterprise™

Notes:

Supports DICOM MWL for order query. ECGs are sent as DICOM waveform objects via a C-STORE transaction.

Website:

www.scimage.com

Other DICOM systems that should be compatible but have not yet been validated:

ChipSoft BV

System Name:

ChipSoft HiX RMR

Cordiana Medical Informatics AG

System Name:

Cordiana Dx16 DICOM ECG Viewer/Editor

Esaote®

System Name:

Suitestensa

NEXUS Enterprise Diagnostics BV (formerly RVC BV)

System Name:

Clinical Assistant

Sectra Medical

System Name:

Sectra PACS

VISUS Technology Transfer GmbH

System Name:

JiveX® ECG

TRADEMARKS

The following marks are the registered trademarks of their respective owners. These marks are used for informational use only and do not imply any endorsement or approval by their respective owners.

<i>Owner</i>	<i>Trademark(s)</i>
Agfa HealthCare, Inc.	<i>HeartStation</i> and <i>Enterprise Imaging</i>
Carestream Health, Inc.	<i>Carestream</i>
Cerner Corporation	<i>PowerChart</i>
Change Healthcare Operations LLC	<i>Change Healthcare</i>
Datamed LLC	<i>Datamed</i>
Edan Instruments, Inc.	<i>Edan</i>
Esaote, S.p.A.	<i>Esaote</i>
Fujifilm Medical Systems USA Inc	<i>Synapse</i>
GE Medical Systems Information Technologies, Inc	<i>MAC, MAC VU360, and MUSE</i>
General Electric Company	<i>GE</i>
Koninklijke Philips Electronics	<i>Philips, PageWriter</i> and <i>IntelliSpace</i>
McKesson Corporation	<i>McKesson</i>
Medimatic S.r.l.	<i>ComPACS</i>
Mortara Instrument, Inc	<i>Mortara, ELI, Burdick, and Pyramis</i>
National Electrical Manufacturers Association Corporation	<i>DICOM</i>
Nihon Kohden Corporation	<i>Cardiofax</i> and <i>Nihon Kohden</i>
Novarad Corporation	<i>NovaCardio</i> and <i>Novarad</i>
Schiller Holding AG	<i>CARDIOVIT</i>
ScImage, Inc	<i>ScImage</i> and <i>PICOM</i>
VidiStar LLC	<i>VidiStar LLC</i>
VISUS Technology Transfer GmbH	<i>JiveX</i>
Welch Allyn, Inc	<i>Welch Allyn</i>