



# IHE Integration Statements

Datamed LLC

DatamedFT™ v2.2

DatamedWL™ v2.2

Doc ID DMD-10000-00200 Rev A

10/2015

### COPYRIGHT

Copyright © 2015 Datamed LLC

Brentwood, TN 37027

(800) 601-3361

+1 901 672 6225

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of Datamed LLC.

### WARRANTY

Datamed LLC ("DLLC") makes no warranty of any kind with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. The software described herein is covered by DLLC's standard EULA, which is available on DLLC's website ([www.datamed.com](http://www.datamed.com)).

### TRADEMARKS

The following marks are the registered trademarks of their respective owners.

<i>Owner</i>	<i>Trademark(s)</i>
Carl E Osborne Jr	<i>Datamed</i>
IHE International, Inc.	<i>IHE</i>

---

## TABLE OF CONTENTS

Table of Contents.....	3
Integration Statement Overview.....	4
DatamedFT™.....	5
DatamedWL™.....	6

---

### INTEGRATION STATEMENT OVERVIEW

IHE Integration Statements describe the conformance of our products with the IHE Technical Framework. This document contains the statements for Datamed® software applications that implement functionality related to the framework.

IHE Integration Statement			
Vendor	Product Name	Version	Date
Datamed LLC	DatamedFT	v2.2	1/1/2014
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Resting ECG Workflow (REWF)	<i>Acquisition Modality</i> when interfacing with a DICOM <i>Integrated ECG Manager</i> . In this configuration, the product is a surrogate for the actual non-DICOM <i>Acquisition Modality</i> and provides an interface for ECG storage. <i>Integrated ECG Manager</i> when interfacing with a DICOM <i>Acquisition Modality</i> . In this configuration, the product is a surrogate for the actual non-DICOM <i>Integrated ECG Manager</i> and provides an interface for ECG storage.	None	
Internet address for Datamed LLC’s IHE information: <a href="http://www.datamed.com/ihe">www.datamed.com/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
DICOM	<a href="http://www.datamed.com/docs/dicom_conformance_statement.pdf">www.datamed.com/docs/dicom_conformance_statement.pdf</a>		
Links to General Information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	

IHE Integration Statement			
Vendor	Product Name	Version	Date
Datamed LLC	DatamedWL	v2.2	6/1/2015
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Resting ECG Workflow (REWF)	<i>Acquisition Modality</i> when interfacing with a DICOM <i>Integrated ECG Manager</i> . In this configuration, the product is a surrogate for the actual non-DICOM <i>Acquisition Modality</i> and provides an interface for the <i>Query Enhanced Modality Worklist</i> transaction.  <i>Integrated ECG Manager</i> when interfacing with a DICOM <i>Acquisition Modality</i> . In this configuration, the product is a surrogate for the actual non-DICOM <i>Integrated ECG Manager</i> and provides an interface for the <i>Query Enhanced Modality Worklist</i> transaction.	None	
Internet address for Datamed LLC’s IHE information: <a href="http://www.datamed.com/ihe">www.datamed.com/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
DICOM	<a href="http://www.datamed.com/docs/dicom_conformance_statement.pdf">www.datamed.com/docs/dicom_conformance_statement.pdf</a>		
Links to General Information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	